

## FLEXURAL TESTING MACHINE

Code : B002



- Used to determine the flexural strength of the test specimen
- Different models are available:
  - C-Frame
    - Load Sensor: Loadcell
  - Rectangular Frame
    - Load Sensor: Pressure Transducer
- Different Capacities are available:
  - 200 kN
  - 300 kN
  - 600 kN
- Can be purchased as stand alone machine or connected to an existing hydraulic unit (must be mentioned at the time of order)

**ORDERING INFORMATION**

Item	Code
FLEXURAL TESTING MACHINE (C-Frame) - 200 kN	B002A20XU
FLEXURAL TESTING MACHINE (C-Frame) - 300 kN	B002A30XU
POWER PACK (for Flexure Frame) - 400 Bar	B002P001U
FLEXURAL C-FRAME - 200 kN	B002P20XH
FLEXURAL C-FRAME - 300 kN	B002P30XH
FLEXURAL TESTING MACHINE (Rectangular Frame) - 200 kN	B003A20XU
FLEXURAL TESTING MACHINE (Rectangular Frame) - 300 kN	B003A30XU
FLEXURAL TESTING MACHINE (Rectangular Frame) - 600 kN	B003A60XU
FLEXURE APPARATUS for CONCRETE BEAMS	B003P001H
FLEXURAL RECTANGULAR FRAME - 200 kN	B003P20XH
FLEXURAL RECTANGULAR FRAME - 300 kN	B003P30XH
FLEXURAL RECTANGULAR FRAME - 600 kN	B003P60XH

**OTHER PHOTOS**

