


## ACCESSORIES for COMPRESSION TESTER

Code : B003



Some of the accessories and spare parts for compression testers series are listed here. Feel free to contact us for further information and guidance.

### **ORDERING INFORMATION**

Item	Code
PISTON SEAL for 1000 kN B-001 and B-004	B001P001H
PISTON SEAL for 1125 kN B-001 and B-004	B001P002H
PISTON SEAL for 1500 kN B-001 and B-004	B001P003H
PISTON SEAL for 2000 kN (for B-001 and B-004)	B001P004H
PISTON SEAL for 3000 kN B-001 and B-004	B001P005H

PISTON SEAL for 5000 kN B-001 and B-004	B001P006H
COMPRESSION FRAME - 1000 kN, Mono-Block, ASTM	B001P007H
COMPRESSION FRAME - 1125 kN, Mono-Block, ASTM	B001P008H
COMPRESSION FRAME - 1500 kN, Mono-Block, ASTM	B001P009H
COMPRESSION FRAME - 2000 kN, Mono-Block, ASTM	B001P010H
COMPRESSION FRAME - 3000 kN, Mono-Block, ASTM	B001P011H
CENTERING APPARATUS for [B 001] SERIES	B001P012H
BELLOWS RUBBER for PRESS SPHERICAL SEATING	B001P013H
COMPRESSION FRAME - 1000 kN, Mono-Block, EN	B001P014H
COMPRESSION FRAME - 1125 kN, Mono-Block, EN	B001P015H
COMPRESSION FRAME - 1500 kN, Mono-Block, EN	B001P016H
COMPRESSION FRAME - 2000 kN, Mono-Block, EN	B001P017H
COMPRESSION FRAME - 3000 kN, Mono-Block, EN	B001P018H
COMPRESSION TESTER UPGRADE PACK	B001P019H
SPHERICAL SEATING for B001	B001P020H
HYDRAULIC VALVE for B001 - 400 bar	B001P021H
SPARE 7" DISPLAY PANEL for B001	B001P022H
ELECTRONIC CONTROLLING UNIT for B001	B001P023H
SPARE LOAD INDICATOR for B001S	B001P024H
SPARE PRESSURE SENSOR (400 bar) for B001	B001P025H
LIMIT SWITCH (with Protective Housing) for B001	B001P026H
SPARE HYDRAULIC PUMP SET (HAWE)	B001P027H
POWER PACK (for Compression Frame) - 400 Bar	B002P000U
POWER PACK (for Flexure Frame) - 400 Bar	B002P001U

POWER PACK (for Compression Frame) - 400 Bar (Without Controlling Unit)	B002P002U
POWER PACK (Servo-Motor, for Compression Frame) - 400 Bar	B002P003U
FLEXURAL C-FRAME - 200 kN	B002P20XH
FLEXURAL C-FRAME - 300 kN	B002P30XH
FLEXURE APPARATUS for CONCRETE BEAMS	B003P001H
RECTANGULAR COMPRESSION PLATES - 30x40 cm	B003P003H
SPACER DISK - 40 mm	B003P004H
SPACER DISK - 50 mm	B003P005H
SPACER DISK - 100 mm	B003P006H
SPLITTING TENSILE APPARATUS for Ø 100 mm CYLINDRER SPECIMENS	B003P007H
SPLITTING TENSILE APPARATUS for Ø 150 mm CYLINDRER SPECIMENS	B003P008H
SPLITTING TENSILE APPARATUS for 150 mm CUBES	B003P009H
THERMAL PRINTER for COMPRESSION TESTER (B-001 and B-004)	B003P010H
METAL DOOR for 2000 kN COMPRESSION TESTER	B003P011H
METAL DOOR for 3000 kN COMPRESSION TESTER	B003P012H
COMPRESSOMETER/EXTENSOMETER (DIAL GAUGE) for Ø75x150mm CONCRETE CYLINDERS	B003P020H
COMPRESSOMETER/EXTENSOMETER (DIGITAL) for Ø75x150mm CONCRETE CYLINDERS	B003P021H
COMPRESSOMETER/EXTENSOMETER (LVDT) for Ø75x150mm CONCRETE CYLINDERS	B003P022H
COMPRESSOMETER/EXTENSOMETER (DIAL GAUGE) for Ø100x200mm CONCRETE CYLINDERS	B003P023H
COMPRESSOMETER/EXTENSOMETER (DIGITAL) for Ø100x200mm CONCRETE CYLINDERS	B003P024H
COMPRESSOMETER/EXTENSOMETER (LVDT) for Ø100x200mm CONCRETE CYLINDERS	B003P025H
COMPRESSOMETER/EXTENSOMETER (DIAL GAUGE) for Ø150x300mm CONCRETE CYLINDERS	B003P026H
COMPRESSOMETER/EXTENSOMETER (DIGITAL) for Ø150x300mm CONCRETE CYLINDERS	B003P027H
COMPRESSOMETER/EXTENSOMETER (LVDT) for Ø150x300mm CONCRETE CYLINDERS	B003P028H

COMPRESSOMETER/EXTENSOMETER (DIAL GAUGE) for Ø160x320mm CONCRETE CYLINDERS	B003P029H
COMPRESSOMETER/EXTENSOMETER (DIGITAL) for Ø160x320mm CONCRETE CYLINDERS	B003P030H
COMPRESSOMETER/EXTENSOMETER (LVDT) for Ø160x320mm CONCRETE CYLINDERS	B003P031H
COMPRESSOMETER/EXTENSOMETER (DIAL GAUGE) for 100x100mm CONCRETE CUBES	B003P032H
COMPRESSOMETER/EXTENSOMETER (DIGITAL) for 100x100mm CONCRETE CUBES	B003P033H
COMPRESSOMETER/EXTENSOMETER (LVDT) for 100x100mm CONCRETE CUBES	B003P034H
COMPRESSOMETER/EXTENSOMETER (DIAL GAUGE) for 150x150mm CONCRETE CUBES	B003P035H
COMPRESSOMETER/EXTENSOMETER (DIGITAL) for 150x150mm CONCRETE CUBES	B003P036H
COMPRESSOMETER/EXTENSOMETER (LVDT) for 150x150mm CONCRETE CUBES	B003P037H
COMPRESSOMETER/EXTENSOMETER (DIAL GAUGE) for 200x200mm CONCRETE CUBES	B003P038H
COMPRESSOMETER/EXTENSOMETER (DIGITAL) for 200x200mm CONCRETE CUBES	B003P039H
COMPRESSOMETER/EXTENSOMETER (LVDT) for 200x200mm CONCRETE CUBES	B003P040H
PISTON for COMPRESSION TESTER (2000 kN)	B003P041H
PISTON for COMPRESSION TESTER (3000 kN)	B003P042H
FLEXURAL RECTANGULAR FRAME - 200 kN	B003P20XH
FLEXURAL RECTANGULAR FRAME - 300 kN	B003P30XH
FLEXURAL RECTANGULAR FRAME - 600 kN	B003P60XH
COMPRESSION FRAME - 1000 kN, 4-Column, EN	B004P001H
COMPRESSION FRAME - 1125 kN, 4-Column, EN	B004P002H
COMPRESSION FRAME - 1500 kN, 4-Column, EN	B004P003H
COMPRESSION FRAME - 2000 kN, 4-Column, EN	B004P004H
COMPRESSION FRAME - 3000 kN, 4-Column, EN	B004P005H
COMPRESSION FRAME - 1000 kN, 4-Column, EN	B004P006H
COMPRESSION FRAME - 1125 kN, 4-Column, EN	B004P007H

---

COMPRESSION FRAME - 1500 kN, 4-Column, EN	B004P008H
COMPRESSION FRAME - 2000 kN, 4-Column, EN	B004P009H
COMPRESSION FRAME - 3000 kN, 4-Column, EN	B004P010H

**OTHER PHOTOS**

