

# CYLINDER MOULD

Code : B010



- Available as:
  - Plastic
  - Steel
- Different dimensions are available

## **ORDERING INFORMATION**

Item	Code
CYLINDER MOULD - Metal (Ø 10 cm)	B010P001H
CYLINDER MOULD - Metal (Ø 15 cm)	B010P002H
CYLINDER MOULD - Metal (Ø 16 cm)	B010P003H

---

CYLINDER MOULD - Metal (Ø 16 cm)	B010P004H
CYLINDER MOULD - Metal (Ø 25 cm)	B010P005H
CYLINDER MOULD - Plastic (Ø 10 cm)	B010P006H
CYLINDER MOULD - Plastic (Ø 15 cm)	B010P007H
CYLINDER MOULD - Plastic (Ø 16 cm)	B010P008H
CYLINDER MOULD (SPRING LID TYPE) - 100 x 200 mm	B010P015H
CYLINDER MOULD (BILLY LID TYPE) - 100 x 200 mm	B010P025H

**OTHER PHOTOS**

