

## CUBE MOULD

Code : B011



- Available as:
  - Plastic
  - Steel
  - Cast Iron
- Different dimensions are available

### ORDERING INFORMATION

Item	Code
CUBE MOULD - Steel (10 cm) : 2-Gang	B011P001H
CUBE MOULD - Cast (15 cm)	B011P002H

CUBE MOULD - Cast (20 cm)	B011P003H
CUBE MOULD - Plastic (10 cm) - 2-Gang	B011P004H
CUBE MOULD - Plastic (15 cm)	B011P005H
COMPACTING BAR (Square Cross-Section) - Steel	B011P006H
COMPACTING BAR (Square Cross-Section) - Wooden	B011P007H

**OTHER PHOTOS**

