


CAPPING PADS & RETAINERS

Code : B016



Unbonded capping offers a safe, precise, and cost-effective alternative to hot sulfur capping compounds for compressive strength testing of concrete specimens.

STANDARDS

ASTM C1231 • AASHTO T22 • AASHTO T851

ORDERING INFORMATION

| Item | Code |
|-----------------------|-----------|
| CAPPING RETAINER - 3" | B016P001H |
| CAPPING RETAINER - 4" | B016P002H |

| | |
|-------------------------------|-----------|
| CAPPING RETAINER - 6" | B016P003H |
| CAPPING PAD - Ø 2" / 70 Shore | B016X27XH |
| CAPPING PAD - Ø 3" / 70 Shore | B016X37XH |
| CAPPING PAD - Ø 4" / 70 Shore | B016X47XH |
| CAPPING PAD - Ø 6" / 70 Shore | B016X67XH |

OTHER PHOTOS

