


## CAPPING PADS & RETAINERS

Code : B016



Unbonded capping offers a safe, precise, and cost-effective alternative to hot sulfur capping compounds for compressive strength testing of concrete specimens.

### **STANDARDS**

ASTM C1231 • AASHTO T22 • AASHTO T851

### **ORDERING INFORMATION**

| Item                  | Code      |
|-----------------------|-----------|
| CAPPING RETAINER - 3" | B016P001H |
| CAPPING RETAINER - 4" | B016P002H |

---

|                               |           |
|-------------------------------|-----------|
| CAPPING RETAINER - 6"         | B016P003H |
| CAPPING PAD - Ø 2" / 70 Shore | B016X27XH |
| CAPPING PAD - Ø 3" / 70 Shore | B016X37XH |
| CAPPING PAD - Ø 4" / 70 Shore | B016X47XH |
| CAPPING PAD - Ø 6" / 70 Shore | B016X67XH |

**OTHER PHOTOS**

