

CONCRETE > SELF-CONSOLIDATING CONCRETE

## V-FUNNEL TEST SET

Code : B025



- The V-Funnel test is a widely used method for assessing the flowability and viscosity of self-compacting concrete (SCC).
- This test helps determine the time it takes for the SCC to flow through a V-shaped funnel, which is indicative of the concrete's ability to pass through tight spaces and around reinforcement without segregation or blockage.

### **STANDARDS**

EN 12350-9

### **TECHNICAL SPECIFICATIONS**

- Made of stainless steel

- Stand-mounted

**ORDERING INFORMATION**

Item	Code
V-FUNNEL TEST SET	B025X00XH