

CONCRETE > SELF-CONSOLIDATING CONCRETE

## **U-SHAPE BOX (SCC)**

Code : B029



- The U-Shaped Box apparatus, specified by EN 11044 standards, is employed to evaluate the flow properties and passing ability of self-compacting concrete (SCC).
- This test assesses how well SCC can flow and fill formwork without segregation, ensuring uniform distribution of aggregates and paste.

### **STANDARDS**

EN 11044 • RILEM REPORT N. 23

### **TECHNICAL SPECIFICATIONS**

- Made of stainless steel
- Guillotine form gate

**ORDERING INFORMATION**

Item	Code
U-SHAPE BOX (SCC)	B029X00XH