

## PULL-OUT TESTER

Code : B106



Used to determine the pullout strength of hardened concrete by measuring the force required to pull an anchor steel rebar from a concrete test specimen or structure. The insert or rebar is either cast into fresh concrete or installed in hardened concrete.

### **SUPPLIED WITH**

- Hydraulic extraction unit (jack and piston assembly)
- Digital display
- Three sets of grippers, suitable to test:
  - Ø 4 mm - Ø 8 mm
  - Ø 10 mm - Ø 20 mm
  - Ø 20 mm - Ø 30 mm

**ORDERING INFORMATION**

Item	Code
PULL-OUT TESTER - 100 kN	B106H10XD
PULL-OUT TESTER - 300 kN	B106H30XD
PULL-OUT TESTER - 500 kN	B106H50XD
GRIPPERS for PULL-OUT TEST - Ø 4-8 mm	B106P001H
GRIPPERS for PULL-OUT TEST - Ø 10-20 mm	B106P002H
GRIPPERS for PULL-OUT TEST - Ø 20-30 mm	B106P003H