

# CONCRETE MIXER

Code : B120



## Pan-Type Mixers

- Used to create trial specimens for the concrete mix design.
- The pan-type mixers provides efficient and homogeneous mixtures.
- 42 lt and 100 lt capacity models are equipped with removable pan that tilts for easier access and emptying on completion of the mixing operation.
- The pan is rotated by a turntable driven by electric motor
- Some models are equipped with secondary motor for double-action mixing (as shown in the table below)

## Drum Type Mixers

- Used to create concrete mixes and specimens.
- The drum-type mixers provides efficient and homogeneous mixtures.
- Equipped with a manual discharge placed at the side of the mixer for direct emptying for the drum.

## **STANDARDS**

EN 1766

## **ORDERING INFORMATION**

Item	Code
PAN TYPE MIXER (Double-acting mixing, 42 / 20 lt)	B120X001U
PAN TYPE MIXER (100 / 56 lt)	B120X002U
PAN for PAN TYPE MIXER (100 / 56 lt)	B120P001H
PAN TYPE MIXER (Double-acting mixing, 100 / 56 lt)	B120X003U
PAN TYPE MIXER (Double-acting mixing, 100 / 56 lt) - 380 V	B120X003V
PAN TYPE MIXER (120 / 100 lt)	B120X004V
PAN TYPE MIXER (285 / 225 lt)	B120X005V
DRUM TYPE MIXER (125 / 100 lt)	B120X006U
DRUM TYPE MIXER (340 / 250 lt)	B120X007V
DRUM TYPE MIXER (490 / 350 lt)	B120X008V
DRUM TYPE MIXER - Diesel Engine (490 / 350 lt)	B120X009H

## **OTHER PHOTOS**

