

VIBRATING TABLE

Code : B125



- Designed from robust steel and has a very steady frame
- Used to compact the concrete specimens in the laboratory
- Equipped with a vibrating motor
- Supplied with a clamp assembly for moulds.

STANDARDS

ASTM C192 • EN 12390-2

ORDERING INFORMATION

| Item | Code |
|-------------------------|-----------|
| VIBRATING TABLE (Small) | B125X001U |

| | |
|--|-----------|
| VIBRATING TABLE (Large) | B125X002U |
| VIBRATING TABLE (Large) - 380V | B125X002V |
| CLAMPING APPARATUS for VIBRATING TABLE (Small) | B125P001U |
| CLAMPING APPARATUS for VIBRATING TABLE (Large) | B125P002U |