

POKER VIBRATOR

Code : B128



- The lab-type poker vibrator is used for internal compaction for the concrete specimens.
- The flexible shaft is 1 meter long.

TECHNICAL SPECIFICATIONS

- Ø 25 mm Flexible shaft (1 m long)
- Direct Connection to the Motor
- 220 - 240 V / 50 - 60 Hz

ORDERING INFORMATION

| Item | Code |
|--------------------------|-----------|
| POKER VIBRATOR - Ø 25 mm | B128X25XU |