

## VEBE CONSISTOMETER

Code : B130



- Used to determine the consistency of concrete using a vibrating table
- Equipped with Vibrating Motor having fixed amplitude and frequency
- ASTM C1176 model is used to prepare cylindrical test specimens from concrete when the standard procedures of rodding and internal vibration are not applicable

### STANDARDS

ASTM C1170 • ASTM 1176 • EN 12350-3

### TECHNICAL SPECIFICATIONS

- Fixed amplitude and frequency
- 220 - 240 V / 50 - 60 Hz

**ORDERING INFORMATION**

Item	Code
VEBE CONSISTOMETER - EN	B130X00EU
VEBE CONSISTOMETER - ASTM	B130X50AU
VEBE VIBRATOR for CYLINDERS - ASTM	B130X20AU