

CONCRETE > SETTING TIME BY PENETRATION

# CONCRETE PENETROMETER

Code : B135



- Used to determine the setting time of concrete, with slump greater than zero, by means of penetration resistance measurements on mortar sieved from the concrete mixture
- Consists of a spring penetrometer with capacity of 100 kgf and interchangeable stainless steel needles
- Equipped with a sliding ring that indicates the load reached upon testing

## **STANDARDS**

ASTM C403

## **TECHNICAL SPECIFICATIONS**

- Stainless Steel Needles Set
- Carrying Case

**ORDERING INFORMATION**

Item	Code
CONCRETE PENETROMETER	B135X00AH
NEEDLES SET for [B 135]	B135P001H