

# CURING TANK

Code : B140



- Used to cure cylinders, cubes, beams and other samples that requires total immersion in water at a specified constant temperature.
- Equipped with a digital thermostat and indicator.
- Supplied with a base rack and a water circulation pump.

## **STANDARDS**

ASTM C31 • BS 1881-11

## **TECHNICAL SPECIFICATIONS**

- 220 - 240 V / 50 - 60 Hz

**SUPPLIED WITH**

- Heater Unit with Digital Thermostat and indicator
- Base Rack
- Water Circulation Pump

**ORDERING INFORMATION**

| Item                                     | Code      |
|--|-----------|
| CURING TANK CIRCULATION PUMP             | B140P001U |
| CURING TANK CONTROLLING UNIT with HEATER | B140P002U |
| BASE RACK 95 x 195 cm for CURING TANK    | B140P003H |
| CURING TANK - Steel                      | B140X001U |
| CURING TANK - Plastic - Large            | B140X002U |
| CURING TANK - Plastic                    | B140X003U |
| CURING TANK - Stainless Steel            | B140X004U |

**OTHER PHOTOS**

