

## **CRACK MICROSCOPE**

Code : B250



- Used to measure the crack width in the concrete elements.
- Supplied with an adjustable light source and adjustable focus to have a better and clearer view for the crack.

### **TECHNICAL SPECIFICATIONS**

- Integral illumination
- Rotating eyepiece
- Focus adjustment
- Battery-Operated

### **ORDERING INFORMATION**

---

Item	Code
CRACK MICROSCOPE (40x)	B250X04XB
CRACK MICROSCOPE (100x)	B250X10XB