

CONCRETE > SETTING TIME BY PENETRATION

# CONCRETE DIAL PENETROMETER

Code : B308



Used to determine the initial set of concrete mortar

## **STANDARDS**

ASTM C403 • AASHTO T197 • UNI 7123

## **TECHNICAL SPECIFICATIONS**

- Capacity: 50 kgf (700 lbs/sq.in = 700 psi)
- Dial with scale readable in both (kgf/cm<sup>2</sup>) and (lbf/sq.in)
- Plunger's surface area: 32.30 mm<sup>2</sup> (~ 1/20 sq.in)

## **ORDERING INFORMATION**

---

| Item                                           | Code      |
|------------------------------------------------|-----------|
| CONCRETE INITIAL SET PENETROMETER - Dial Model | B308X00XH |