


## VIDEO EXTENSOMETER

Code : E023



- Ability to analyze 4 points on the specimen and measure the axial strain between the given points
- High accuracy confirming to the international standards
- Does not need to be clipped, glued or connected to the specimen
- Wide range of specimen (cubes, cylinder, rods, flat, strips, irregular shpes ... etc)
- Used to measure: strain, deformation, elongation, deflection, modulus of elasticity ... etc combined with ALFA Equipment and ALFA State-of-the-art software series

### **SUPPLIED WITH**

- Tripod Stand
- LED Light
- Marker Pen
- Marking Stickers

- Connection Cables
- PC Software

**ORDERING INFORMATION**

| Item               | Code      |
|--------------------|-----------|
| VIDEO EXTENSOMETER | E023X001X |