

## CRADLE

Code : G100P001H



- Used to hang the specimens under weighing balance in Specific Gravity Frame to determine the mass of specimen in the air and inside the water to calculate the specific gravity of the specimens
- Made of stainless-steel rod
- To be used with Specific Gravity Frame

### **TECHNICAL SPECIFICATIONS**

- Made from Stainless-steel

### **ORDERING INFORMATION**

Code	Name
G 100 X00X H	SPECIFIC GRAVITY FRAME

G 100 P001 H

CRADLE

G 100 P002 H

DENSITY BASKET (Wire Basket)

G 100 P003 H

WATER TANK

G 100 P004 H

SOUNDNESS BASKET (Mesh Size 1 mm)