

VACUUM PUMP

Code : G702



Used to vacuum specimens, remove air voids and speed up deairing process.

AVAILABLE MODELS

Code	Item
G 702 X051 U	VACUUM PUMP - 51 lt/min
G 702 X100 U	VACUUM PUMP - 100 lt/min
G 702 X128 U	VACUUM PUMP - 128 lt/min

ORDERING INFORMATION

Code	Item
G 702 X051 U	VACUUM PUMP - 51 lt/min

G 702 X100 U VACUUM PUMP - 100 lt/min

G 702 X128 U VACUUM PUMP - 128 lt/min

G 702 P001 H HOSE for VACUUM PUMP