


## VERNIER CALIPER

Code : G709



- Used to measure the dimensions of an object, generally by placing two movable points of the instrument across the object or span to be measured.
- Many types of calipers permit reading out a measurement on a ruled scale, a dial, or an electronic digital display. A common association is to calipers using a sliding vernier scale.

### **ORDERING INFORMATION**

Code	Name
G 709 X001 B	DIGITAL CALIPER - 150 mm
G 709 X002 B	DIGITAL CALIPER - 200 mm
G 709 X003 B	DIGITAL CALIPER - 300 mm

---

G 709 X004 B                      DIGITAL CALIPER - 500 mm

G 709 X005 H                      VERNIER CALIPER - 150 mm

G 709 X006 H                      VERNIER CALIPER - 300 mm

G 709 X007 H                      VERNIER CALIPER - 500 mm

**OTHER PHOTOS**

