

## DIAL INDICATOR

Code : G940



- Dial gauges serve a crucial purpose in the realm of measurements, as they are adept at quantifying not only vertical but also horizontal deformations. Their utility spans across a diverse array of tests, encompassing a wide spectrum of industries and applications. By facilitating the precise assessment of both vertical and horizontal shifts or distortions, these gauges contribute indispensably to ensuring accurate data collection and informed decision-making within various testing scenarios.

### ORDERING INFORMATION

Code	Name
G 940 P001 B	DIGITAL DIAL GAUGE - 12.7 / 0.001 mm
G 940 P002 B	DIGITAL DIAL GAUGE - 25 / 0.01 mm

---

G 940 P003 H	DIAL GAUGE - 10 / 0.01 mm
G 940 P004 H	DIAL GAUGE - 12.7 / 0.002 mm
G 940 P005 H	DIAL GAUGE - 25 / 0.01 mm
G 940 P006 H	DIAL GAUGE - 30 / 0.01 mm
G 940 P007 H	DIAL GAUGE - 50 / 0.01 mm

**OTHER PHOTOS**

