

# CALIBRATION WEIGHT

Code : G950



- These weight measurements play a crucial role in conducting regular and periodic checks of the balances. By utilizing these weights, the process of balance verification and validation is not only facilitated but also enhanced. These weight values serve as fundamental reference points against which the existing balances are compared during routine examinations.

## **ORDERING INFORMATION**

Code	Name
G 951 X001 H	CALIBRATION WEIGHT SET - 2110 g
G 950 X001 H	CALIBRATION WEIGHT - 1 g
G 950 X005 H	CALIBRATION WEIGHT - 5 g
G 950 X010 H	CALIBRATION WEIGHT - 10 g

---

G 950 X01K H	CALIBRATION WEIGHT - 1 kg
G 950 X020 H	CALIBRATION WEIGHT - 20 g
G 950 X02K H	CALIBRATION WEIGHT - 2 kg
G 950 X050 H	CALIBRATION WEIGHT - 50 g
G 950 X05K H	CALIBRATION WEIGHT - 5 kg
G 950 X100 H	CALIBRATION WEIGHT - 100 g
G 950 X10K H	CALIBRATION WEIGHT - 10 kg
G 950 X20K H	CALIBRATION WEIGHT - 20 kg
G 950 X25K H	CALIBRATION WEIGHT - 25 kg
G 950 X500 H	CALIBRATION WEIGHT - 500 g