


## JAW CRUSHER

Code : M412



- Used to crush samples with medium-hard, hard, brittle and tough materials. Efficiency and safety makes it ideal for sample preparation in laboratories and industrial plants.
- Application examples are so wide such as construction materials, glass, granite, minerals, quartz, rocks, silicon, slag, coal, oxide ceramics ... etc.
- The materials after crushing have high degree of size reduction and high fineness.
- Double action of jaws that opens/closes from the bottom and the crushing mechanism is from top to bottom.
- Supplied complete with electrical motor.

### AVAILABLE MODELS

Code	Item
M 412 X001 U	JAW CRUSHER

M 412 X002 V                      JAW CRUSHER (Large Capacity)

M 412 X003 V                      JAW CRUSHER (Extra Large Capacity)

**ORDERING INFORMATION**

Code	Item
M 412 X001 U	JAW CRUSHER
M 412 X002 V	JAW CRUSHER (Large Capacity)
M 412 X003 V	JAW CRUSHER (Extra Large Capacity)
M 412 P002 H	SPARE SIDE PLATES for G-412

**OTHER PHOTOS**

