

## **AGGREGATE SHAPE GAUGE**

Code : R036



- Used to measure the shape index, ratio between the length to the thickness of any individual aggregate particle.
- Individual particles in a sample of coarse aggregate are classified based on the ratio of their length to thickness.
- The shape index is calculated as the mass of particles with a ratio of dimensions Length/Thickness more than 3 expressed as a percentage of the total dry mass of particles tested.

### **STANDARDS**

EN 933-4 • EN 933-5 • EN 933-7

### **TECHNICAL SPECIFICATIONS**

Made of steel

**ORDERING INFORMATION**

| Code         | Item                      |
|--------------|---------------------------|
| R 036 X00E H | SHAPE GAUGE (Shape Index) |