

## SOIL PENETROMETER (Pocket Type)

Code : T304



Used to determine the classification of cohesive soil based on its consistency, shear strength and the approximate unconfined compressive strength

### TECHNICAL SPECIFICATIONS

Capacity	Surface Area (mm <sub>c</sub> )	Surface Area (cm <sub>c</sub> )
4.5 kgf/cm <sub>c</sub>	32.30	0.3230
16 kgf/cm <sub>c</sub>	32.30	0.3230

### ORDERING INFORMATION

Item	Code
POCKET PENETROMETER - 4.5 kgf/cm <sub>c</sub>	T 304 X05X H
POCKET PENETROMETER - 16 kgf/cm <sub>c</sub>	T 304 X16X H