

SOIL PENETROMETER (Pocket Type)

Code : T304



Used to determine the classification of cohesive soil based on its consistency, shear strength and the approximate unconfined compressive strength

TECHNICAL SPECIFICATIONS

Capacity	Surface Area (mm _c)	Surface Area (cm _c)
4.5 kgf/cm _c	32.30	0.3230
16 kgf/cm _c	32.30	0.3230

ORDERING INFORMATION

Item	Code
POCKET PENETROMETER - 4.5 kgf/cm _c	T 304 X05X H
POCKET PENETROMETER - 16 kgf/cm _c	T 304 X16X H