

SOIL AUGER SET

Code : T305



- Used to take samples from the surface or near-surface soils. It is also used to extract disturbed samples from a soil bulk for boring logs
- Consists of a T handle, extension tube and an auger head
- The auger head comes in different diameters, to be chosen depending on the soil type and the desired amount of soil to be extracted.

STANDARDS

ASTM D420 • ASTM D1452 • AASHTO T86 • AASHTO T202

ORDERING INFORMATION

| Item | Code |
|-------------------------|--------------|
| SOIL AUGER SET - Ø 2.5" | T 305 X063 H |

| | |
|-----------------------------------|--------------|
| SOIL AUGER SET - Ø 80 mm | T 305 X080 H |
| SOIL AUGER SET - Ø 100 mm | T 305 X100 H |
| SOIL AUGER SET - Ø 150 mm | T 305 X150 H |
| T-HANDLE for AUGER SET | T 305 P005 H |
| EXTENSION ROD for AUGER SET - 1 m | T 305 P006 H |
| AUGER HEAD - Ø 2.5" | T 305 P001 H |
| AUGER HEAD - Ø 80 mm | T 305 P002 H |
| AUGER HEAD - Ø 100 mm | T 305 P003 H |
| AUGER HEAD - Ø 150 mm | T 305 P004 H |