

## SOIL PENETROMETER (Dial Type)

Code : T308



- Used to determine the angle of friction, the cohesion index and the approximate unconfined compressive strength.
- Can be used on the foundation soil from clayey soil to sandy soil.
- The Zero Setting for the equipment can be done by pushing the reset button.

### **SUPPLIED WITH**

Plunger Diameter	Surface Area (mm <sup>2</sup> )	Surface Area (cm <sup>2</sup> )
Ø 6.35 mm	31.67	0.3167
Ø 10 mm	78.54	0.7854
Ø 15 mm	176.71	1.7671
Ø 20 mm	314.16	3.1416

Ø 25 mm	490.87	4.9087
---------	--------	--------

**ORDERING INFORMATION**

Item	Code
SOIL DIAL PENETROMETER SET	T 308 X00X H
PLUNGER for T308 - Ø 6.35 mm	T 308 P001 H
PLUNGER for T308 - Ø 10 mm	T 308 P002 H
PLUNGER for T308 - Ø 15 mm	T 308 P003 H
PLUNGER for T308 - Ø 20 mm	T 308 P004 H
PLUNGER for T308 - Ø 25 mm	T 308 P005 H