


# MARSHALL COMPACTOR, Manual

Code : A016



- Used to determine the resistance of an asphalt specimen to plastic flow.
- The test is performed using a manually dropped hammer.

## **STANDARDS**

ASTM D6926

## **TECHNICAL SPECIFICATIONS**

Composed of:

- Compaction Hammer
- Wooden Pedestal
- Steel Plate

- Mould Holder
- Hammer Guide

**ORDERING INFORMATION**

Item	Code
MARSHALL COMPACTOR, HAND-OPERATED - ASTM	A016H04AH

**OTHER PHOTOS**

