


DISTILLATION for CUTBACK ASPHALT

Code : A450



- Used to determine the amount of water in bitumen and petroleum products. The mechanism of the test is done via distilling the specimen with a water immiscible (volatile solvent).

STANDARDS

ASTM D402 • AASHTO T-78

SET

- Distillation Flask
- Condenser
- Ring Burner
- Thermometer (2ea)

ORDERING INFORMATION

| Item | Code |
|--|-----------|
| DISTILLATION UNIT for CUTBACK ASPHALT - ASTM | A450X00AH |