

WATER in BITUMINOUS MATERIALS by DISTILLATION TEST SET

Code : A455



Used to determine the water content in the range from 0 to 25 % volume in petroleum products, tars, and other bituminous material

STANDARDS

ASTM D95 • ASTM D244 • AASHTO T55

SUPPLIED WITH

- Glass balloon - 500 ml
- Glass receiver - 10 ml
- Glass reflux condenser

- Electrical heater with clamping system

ORDERING INFORMATION

| Item | Code |
|--|-----------|
| WATER in BITUMINOUS MATERIALS by DISTILLATION TEST SET | A455X00AH |