

# FLASH and FIRE POINT APPARATUS (Cleveland Open Cup)

Code : A470



- Used to determine the maximum temperature bitumen can be safely heated without the danger of an instantaneous flash in the presence of an open flame.
- Supplied with an electronic heater, a thermoregulator, two thermometers, and a sample cup.

## **STANDARDS**

ASTM D117 • ASTM D92 • AASHTO T-48

## **SUPPLIED WITH**

- Electric heater

- Thermoregulator
- Two Thermometers
- Sample Cup

**ORDERING INFORMATION**

Item	Code
FLASH and FIRE POINT APPARATUS (Cleveland Open Cup)	A470H00XU
AUTOMATIC FLASH and FIRE POINT APPARATUS (Cleveland Open Cup)	A470A00XU