

## CORING MACHINE, Electrical

Code : B061



- Used to take core samples from irregular rocks or to extrude the core from specimens for test purposes.
- Standard 1 ¼" thread
- Core Bits should be ordered separately

### **TECHNICAL SPECIFICATIONS**

- Power input : 2500 W
- Rated Speed : 450 rpm
- Max. drill : Ø 162 mm
- Water inlet
- B-061/R model can rotate from 0° to 89°
- Standard 1 ¼" thread for bits
- 220 - 240 V / 50 - 60 Hz

**ORDERING INFORMATION**

Item	Code
CORING MACHINE	B061M01XU
ROTATABLE CORING MACHINE	B061M02XU
CORE BIT for CONCRETE - $\varnothing$ 50 mm x 400 mm	B061P001H
CORE BIT for CONCRETE - $\varnothing$ 75 mm x 400 mm	B061P002H
CORE BIT for CONCRETE - $\varnothing$ 100 mm x 400 mm	B061P003H
CORE BIT for CONCRETE - $\varnothing$ 150 mm x 400 mm	B061P004H

**OTHER PHOTOS**

