

CONCRETE > SETTING TIME BY PENETRATION

# CONCRETE POCKET PENETROMETER

Code : B304



- Used to determine the initial set of concrete mortar

## **STANDARDS**

ASTM C403 • AASHTO T197 • UNI 7123

## **TECHNICAL SPECIFICATIONS**

- Capacity: 4.5 MPa
- Plunger's surface area: 32.30 mm<sup>2</sup> (~ 1/20 sq.in)

## **ORDERING INFORMATION**

---

Item	Code
CONCRETE POCKET PENETROMETER - Pen-Type	B304X00XH