


ACCESSORIES for CEMENT TESTER

Code : C003



- To be used with Cement Compression/Flexure machines.
- Jigs are used to apply homogenous load on specimens as per the international standards

STANDARDS

ASTM C109 • ASTM C349 • EN 196-1

ORDERING INFORMATION

| Item | Code |
|--|-----------|
| DUCTILITY TESTER with COOLING UNIT [60 Hz] | A440C00XK |
| DUCTILITY TESTER with COOLING UNIT | A440C00XT |

| | |
|--|-----------|
| CEMENT COMPRESSION and FLEXURE TESTER - TFT | C001A25XU |
| HIGH-CAPACITY CEMENT COMPRESSION and FLEXURE TESTER - TFT | C001A60XU |
| CEMENT COMPRESSION TESTER - 250 kN | C002A25XU |
| CEMENT COMPRESSION TESTER - 600 kN | C002A60XU |
| COMPRESSION JIG - 40 mm | C003P001H |
| COMPRESSION JIG - 50 mm | C003P002H |
| FLEXURE JIG | C003P003H |
| COMPRESSION JIG - 70.7 mm | C003P004H |
| TRANSVERSE DEFORMATION TESTER - 250 N | C004A250U |
| HIGH-CAPACITY CEMENT COMPRESSION and FLEXURE TESTER (Loadcell Model) | C005A60XU |

OTHER PHOTOS

