


LENGTH COMPARATOR

Code : C018



- Used to determine the length change in hardened cement paste, mortar, and concrete prisms
- Equipped with Digital Gauge - 12.7 / 0.001 mm (DI/D-12/0.001)
- Height adjustable to suit specimens up-to 300 mm

STANDARDS

ASTM C151 • ASTM C490 • EN 12617-4 • EN 1367-4

ORDERING INFORMATION

Item	Code
DIGITAL LENGTH COMPARATOR - Digital	C018X01XB
ANALOG LENGTH COMPARATOR - Analog	C018X02XH

REFERENCE ROD - 160 mm	C018P001H
REFERENCE ROD - 200 mm	C018P002H
REFERENCE ROD - 250 mm	C018P003H
REFERENCE ROD - 254 mm	C018P004H
REFERENCE ROD - 280 mm	C018P005H