


# LE CHATELIER SOUNDNESS TEST SET

Code : C025



- Used to determine the soundness of the cement.

## **STANDARDS**

EN 196-3

## **SET**

- Le chatelier moulds (3 ea)
- Glass plates - 50 x 50 mm (6 ea)
- Weight - 100 g (3 ea)
- Tamping Rod - Ø 17 mm x 70 g (1 ea)
- Steel Ruler (1 ea)

**ORDERING INFORMATION**

Item	Code
LE CHATELIER MOULD	C023X00EH
LE CHATELIER SOUNDNESS TEST SET	C025X00EH
LE CHATELIER BATH - 6 lt	G040X06XU
LE CHATELIER FLASK - 250 ml	W019G250H
LE CHATELIER STAND for [G 040 X06X U]	G040P001H
SPARE HEATER for LE CHATELIER BATH	G040P004H