

CEMENT MIXER

Code : C050



- Used to assure an efficient mixing for the cement and the mortar samples.
- Can take a capacity up to 5 liters.
- Equipped with a planetary mixing mechanism for better performance.
- The beater speed is adjustable between 140 rpm and 285 rpm.

STANDARDS

EN 196-1 • EN 196-3

TECHNICAL SPECIFICATIONS

- Automatic mixing

- Bowl Capacity: 5 Lt.
- Planetary mixing action.
- Beater Speeds: 140 rpm / 285 rpm.
- 220 - 240 V / 50 - 60 Hz

ORDERING INFORMATION

Item	Code
AUTOMATIC PROGRAMMABLE CEMENT MIXER	C050A00XU
MORTAR MIXER	C050M00XU
BEATER for MORTAR MIXER	C050P001H
STEEL BOWL for MORTAR MIXER	C050P002H
BEATER (Stainless Steel) for MORTAR MIXER	C050P003H

OTHER PHOTOS

