


# CEMENT VIBRATING MACHINE

Code : C126



- Used to vibrate cement mortar cube sample in the 70.7 mm cube moulds.
- Mould should be ordered separately

## **STANDARDS**

BS 4550

## **TECHNICAL SPECIFICATIONS**

- Vibrates at 12,000 cycles per minute
- 220 - 240 V / 50 - 60 Hz

**ORDERING INFORMATION**

Item	Code
VIBRATING MACHINE for 70.7 mm CUBE MOULDS	C126X00BT
VIBRATING MACHINE for 70.7 mm CUBE MOULDS [60 Hz]	C126X00BK