

## LARGE CAPACITY SAMPLE SPLITTER (Universal)

Code : G081



- Designed to halve/divide aggregates, soils, sands, and gravel that are too large in volume into two representative halves.
- The lever-controlled unit is made of heavy-duty steel with a hopper holding up to 28 liters sample.
- The sample is divided into two equal parts using a set of chutes. These chutes release the material in alternating directions into individual collecting pans.
- Each adjustable chute having a thickness of 12.5 mm provides flexibility, allowing to have 12.5 mm to 150 mm opening sizes.

### **STANDARDS**

ASTM C136 • EN 933-3 • AASHTO T27

### **TECHNICAL SPECIFICATIONS**

- Made of heavy-duty steel
- Coating for rust protection
- Adjustable Chutes
  - Thickness: 12.5 mm
  - Count: 48
- Hopper Capacity : 28 liters
- Collecting Pan : 22 liters
- Adjustable openings from 12.5 mm to 150 mm at 12.5 mm intervals

### **SUPPLIED WITH**

2 Collecting pans (Each 22 liter)

### **ORDERING INFORMATION**

Code	Item
G 081 X00X H	UNIVERSAL SAMPLE SPLITTER
G 081 P001 H	COLLECTING PAN for UNIVERSAL SAMPLE SPLITTER