


# VERNIER CALIPER

Code : G709



- Used to measure the dimensions of an object, generally by placing two movable points of the instrument across the object or span to be measured.
- Many types of calipers permit reading out a measurement on a ruled scale, a dial, or an electronic digital display. A common association is to calipers using a sliding vernier scale.

## **ORDERING INFORMATION**

Item	Code
DIGITAL CALIPER - 150 mm	G709X001B
DIGITAL CALIPER - 200 mm	G709X002B
DIGITAL CALIPER - 300 mm	G709X003B

---

DIGITAL CALIPER - 500 mm	G709X004B
VERNIER CALIPER - 150 mm	G709X005H
VERNIER CALIPER - 300 mm	G709X006H
VERNIER CALIPER - 500 mm	G709X007H
DIGITAL MICROMETER - 25 mm	G709X101B

**OTHER PHOTOS**

