

# HOTPLATE

Code : G711



- A hotplate in engineering testing laboratories is used to heat samples or solutions uniformly during tests. It is essential for procedures like drying, mixing, or chemical reactions, ensuring controlled temperatures and consistent results.

## **ORDERING INFORMATION**

Item	Code
RECTANGULAR HOTPLATE - 30 x 30 cm	G711X001U
RECTANGULAR HOTPLATE - 35 x 45 cm	G711X002U
RECTANGULAR HOTPLATE - 40 x 60 cm	G711X003U
HOTPLATE w/ MAGNETIC STIRRER	G711X004U

---

ROUND HOTPLATE - Ø 170 mm

G711X005U

ROUND HOTPLATE - Ø 170 mm (2-Place)

G711X006U

**OTHER PHOTOS**

