


## WIRE BRUSH

Code : G718



- A wire brush in engineering testing laboratories is also used to clean sieves by removing trapped particles from the mesh. This helps maintain sieve accuracy, ensures proper gradation results, and prolongs the lifespan of the sieves.

### **ORDERING INFORMATION**

Item	Code
WIRE BRUSH	G718X001H
HAIR BRUSH - No: 1	G719X001H
HAIR BRUSH - No: 2	G719X002H
HAIR BRUSH - No: 3	G719X003H
HAIR BRUSH - No: 4	G719X004H

BRUSH with HANDLE

G719X005H