

ARMORED THERMOMETER w/ CASE

Code : G781



- Thermometers in engineering testing laboratories are used to measure the temperature of samples, materials, or equipment during experiments. They help ensure accurate temperature control and monitoring, which is critical for processes that are temperature-sensitive, such as curing, heating, or material testing.

ORDERING INFORMATION

Item	Code
GLASS THERMOMETER - 100°C	G780X100H
GLASS THERMOMETER - 200°C	G780X200H
GLASS THERMOMETER - 250°C	G780X250H
GLASS THERMOMETER - 360°C	G780X360H

ARMORED THERMOMETER w/ CASE - 200°C	G781X200H
DIGITAL THERMOMETER - 300°C	G782X300B
THERMO-HYGRO METER	G783X001B
DIAL THERMOMETER (-40°C / +70°C)	G784X070H
DIAL THERMOMETER (0°C / +200°C)	G784X200H
INFRARED THERMOMETER	G785X001B
MIN-MAX THERMOMETER (-40°C / +50°C)	G786X001B