

NUCLEAR DENSITY/MOISTURE TESTER

Code : G852



- Used to determine density and the moisture content of the sample in a very short period with very accurate results.
- The testing time can be chosen from the control panel as 15 seconds, 1 minute and 4 minutes.
- The depth range that can be measured by this equipment is 12" (30 cm), with 1" measurement increment.

TECHNICAL SPECIFICATIONS

- Test time: 15 sec / 1 min / 4 min
- Depth Range up to 12" (30 cm)
- Increment: 1" (2.5 cm)

ORDERING INFORMATION

Item	Code
NUCLEAR DENSITY/MOISTURE TESTER	G852X00AB