

LOADCELL

Code : G905



- Used for load measurements and calibration purposes
- Conforms to ASTM and EN Class 1
- Suitable for both laboratory and field use
- Supplied with durable carrying case

SUPPLIED WITH

- Touch Screen Digital Indicator:
 - Unit Selection (kN / N / kgf / lbf)
 - Taring / Zeroing Option
 - Calibration (Password-Protected)
 - Data can be transferred to PC for reporting and printing in real-time
 - Battery and Mains operated

- Durable Carrying Case
- Calibration Certificate (Internationally traceable)

ORDERING INFORMATION

Code	Name
G 905 X005 D	LOADCELL with DIGITAL INDICATOR - 5 kN
G 905 X01k D	LOADCELL with DIGITAL INDICATOR - 1000 kN
G 905 X025 D	LOADCELL with DIGITAL INDICATOR - 25 kN
G 905 X02K D	LOADCELL with DIGITAL INDICATOR - 2000 kN
G 905 X03K D	LOADCELL with DIGITAL INDICATOR - 3000 kN
G 905 X050 D	LOADCELL with DIGITAL INDICATOR - 50 kN
G 905 X250 D	LOADCELL with DIGITAL INDICATOR - 250 kN
G 905 X500 D	LOADCELL with DIGITAL INDICATOR - 500 kN

OTHER PHOTOS

