

## LOAD-RING (Proving Ring)

Code : G906



- Used for load measurements and calibration purposes
- Suitable for both laboratory and field use
- Supplied with factory calibration report and durable carrying case

### SUPPLIED WITH

- Factory Calibration Report
- Durable Carrying Case

### ORDERING INFORMATION

Code	Name
G 906 X005 H	LOAD RING - 5 kN

---

G 906 X010 H	LOAD RING - 10 kN
G 906 X030 H	LOAD RING - 30 kN
G 906 X050 H	LOAD RING - 50 kN