

# ALUMINUM SAMPLING MACHINE

Code : M062



- Used to prepare circular samples from aluminum ingots for XRF analysis
- Fully automatic, with 3 functions:
  - Drilling
  - Polishing
  - Cutting
- Ingot is placed on the support system and all the a/m functions are done automatically in a short time by an operator
- The sample taken has a desired diameter and thickness with smooth surface

## **TECHNICAL SPECIFICATIONS**

- Fully automatic, with 3 functions:
  - Drilling

- Polishing
- Cutting
- Ingot is placed on the support system and all the a/m functions are done automatically in a short time by an operator

**ORDERING INFORMATION**

Item	Code
ALUMINUM SAMPLING MACHINE	M062X00XV