

## DISC MILL

Code : M420



- Used under rough conditions for grinding of hard, medium-hard, brittle and tough materials
- Required few minutes to achieve the desired grind size, and work continuously
- Application examples: minerals, ores, bauxite, clinkers, chalk, chamotte, concrete, construction materials, soil samples, gypsum ... etc.

### **TECHNICAL SPECIFICATIONS**

Small Disc Model:

- Capacity: 20 – 100 kg/hr depending on material and final size
- Disc diameter: 170 mm
- Feeding size: < 10 mm
- Final fineness: 50 – 200 microns

- Collector: 25 x 20 x 8 cm
- 380 V / 50 Hz

Large Disc Model:

- Capacity: 40 – 200 kg/hr depending on material and final size
- Disc diameter: 205 mm
- Feeding size: < 15 mm
- Final fineness: 50 – 200 microns
- Collector: 25 x 20 x 8 cm
- 380 V / 50 Hz

### **ORDERING INFORMATION**

Item	Code
DISC MILL - Large	M420X001U
DISC MILL - Small	M420X002U
SPARE DISCS for G-420/L	M420P001H
SPARE DISCS for G-420/S	M420P002H

### **OTHER PHOTOS**

